


BestSling User Guide

STAND ASSIST SLINGS




WARNING!

- Step 1: Check sling weight capacity for the patient being lifted. Do not exceed maximum safe working load.
- Step 2: Review following procedures carefully, prior to attempting lift. Call with any questions.
- Step 3: Position sling behind seated patient, with the 2 red pads under arms (Diagram 2).
- Step 4: Cross padded belt over chest/stomach area and secure with snap buckle belt (Diagram 2).
- Step 5: Loop strap B through black strap C and hook to corresponding A & B hooks on mast (Diagram 1 and 4).
- Step 6: Have patient grasp handle bars and place feet on lift platform, with knees against black knee pad (Diagram 5).
- Step 7: Keep arms outside all loops, and lean head slightly back (Diagram 5).
- Step 8: Do not lock caster brakes to allow lift to properly balance weight.
- Step 9: Begin the lift and transfer as in Diagram 6.

* Optional buttock strap shown in Diagram 6 offers extra support and attaches to extra loops on spreader bar below a and b.

SLING LOOP OPTIONS

- Long Loop- reclined position.
- Center Loop- semi-reclined position
- Short Loop- most upright position



1-877-822-9033

