









Spreader Bar Compatibility Chart

Model	Part No.	Lift Model Used On	Lift Model Compatible With	Note
 <p>2 pt. Spreader Bar (Red)</p>	WP-SBAR-PL2RW	PL400EF	PL350CT PL182 PL228 PL273	It can be used on all PL lifts if don't mind the color difference. And will have different lifting range on all LEGACY lifts
 <p>2 pt. Spreader Bar (Red)</p>	WP-PL-228-ML	PL228-ML	PL350CT PL400EF PL182 PL273	It can be used on all PL lifts if don't mind the color difference. And will have different lifting range on all LEGACY lifts
 <p>6 pt. Spreader Bar (Blue)</p>	WP-SBAR-PL6BG	PL400 PL500 PL600		It can be used on all PL lifts if don't mind the color difference. And will have different lifting range on all LEGACY lifts
 <p>6 pt. Spreader Bar (Red)</p>	WP-SBAR-PL6RW	Build for PL400EF	PL350CT PL182 PL228 PL273	It can be used on all PL lifts if don't mind the color difference. And will have different lifting range on all LEGACY lifts

Spreader Bar Compatibility Chart

 <p style="text-align: center;">6 pt. Spreader Bar</p>	<p style="text-align: center;">WP-PL228-011</p>	<p style="text-align: center;">PL182 PL228 PL273</p>	<p style="text-align: center;">PL350CT PL400EF</p>	<p>It can be used on all PL lifts if don't mind the color difference.</p> <p>And will have different lifting range on all LEGACY lifts</p>
	<p style="text-align: center;">WP-PL400H-SBAR</p>	<p style="text-align: center;">PL400H PL400HE</p>	<p style="text-align: center;">PL350CT PL400EF PL182</p>	<p style="text-align: center;">Designed for 400lb ONLY</p>
	<p style="text-align: center;">WP-PL350CT-UBAR</p>	<p style="text-align: center;">PL350CT</p>		<p>Only compatible with Cradle Clip slings</p>
 <p style="text-align: center;">Non-rotating 6 pt. Spreader Bar</p>	<p style="text-align: center;">WP-SBAR-PL400</p>	<p style="text-align: center;">PL400E</p>		<p>It can be used on all PL lifts if don't mind the color difference.</p> <p>And will have different lifting range on all ProCare and Care patien lifts</p>

Spreader Bar Compatibility Chart

 <p style="text-align: center;">Non-rotating 6 pt. Spreader Bar</p>	<p style="text-align: center;">WP-SBAR-PL600</p>	<p style="text-align: center;">PL500E PL600E</p>	<p style="text-align: center;">PL400E</p>	<p>It can be used on all PL lifts if don't mind the color difference.</p> <p>And will have different lifting range on all ProCare and Care patient lifts</p>
 <p style="text-align: center;">2 Pt. Spreader Bar</p>	<p style="text-align: center;">DISCONTINUED</p>	<p style="text-align: center;">PL400E PL500E PL600E</p>		<p>Replace with WP-SBAR-PL400 or WP-SBAR-PL600</p>
 <p style="text-align: center;">Rotating 6 pt. Spreader Bar</p>	<p style="text-align: center;">DISCONTINUED</p>	<p style="text-align: center;">PL400E PL500E PL600E</p>		<p>Replace with WP-SBAR-PL400 or WP-SBAR-PL600</p>
 <p style="text-align: center;">Old Rotating 6 pt. Spreader Bar</p>	<p style="text-align: center;">DISCONTINUED</p>	<p style="text-align: center;">PL400E PL500E PL600E</p>		<p>Replace with WP-SBAR-PL400 or WP-SBAR-PL600</p>