












Spreader Bar Compatibility Chart

Model	Part No.	Lift Model Used On	Lift Model Compatible With	Note
 <p data-bbox="201 570 489 602">2 pt. Spreader Bar (Red)</p>	WP-SBAR-PL2RW	PL400EF	PL350CT PL182 PL228	It can be used on all PL lifts if don't mind the color difference. And will have different lifting range on all LEGACY lifts It is compatible with PL350CT but won't function properly.
 <p data-bbox="201 886 489 919">2 pt. Spreader Bar (Red)</p>	WP-PL-228-ML	PL228-ML	PL350CT PL400EF PL182	It can be used on all PL lifts if don't mind the color difference. And will have different lifting range on all LEGACY lifts It is compatible with PL350CT but won't function properly.
 <p data-bbox="201 1192 489 1224">6 pt. Spreader Bar (Blue)</p>	WP-SBAR-PL6BG	PL400 PL500 PL600		It can be used on all PL lifts if don't mind the color difference. And will have different lifting range on all LEGACY lifts It is compatible with PL350CT but won't function properly.

Spreader Bar Compatibility Chart

 <p style="text-align: center;">6 pt. Spreader Bar (Red)</p>	<p style="text-align: center;">WP-SBAR-PL6RW</p>	<p style="text-align: center;">Build for PL400EF</p>	<p style="text-align: center;">PL350CT PL182 PL228</p>	<p>It can be used on all PL lifts if don't mind the color difference.</p> <p>And will have different lifting range on all LEGACY lifts</p> <p>It is compatible with PL350CT but won't function properly.</p>
 <p style="text-align: center;">6 pt. Spreader Bar</p>	<p style="text-align: center;">WP-PL228-011</p>	<p style="text-align: center;">PL182 PL228 PL273</p>	<p style="text-align: center;">PL350CT PL400EF</p>	<p>It can be used on all PL lifts if don't mind the color difference.</p> <p>And will have different lifting range on all LEGACY lifts</p> <p>It is compatible with PL350CT but won't function properly.</p>
	<p style="text-align: center;">WP-PL400H-SBAR</p>	<p style="text-align: center;">PL400H PL400HE</p>	<p style="text-align: center;">PL400EF PL182</p>	<p style="text-align: center;">Designed for 400lb ONLY</p>
	<p style="text-align: center;">WP-PL350CT-UBAR</p>	<p style="text-align: center;">PL350CT</p>		<p style="text-align: center;">Only compatible with Cradle Clip slings</p>

Spreader Bar Compatibility Chart

 <p style="text-align: center;">Non-rotating 6 pt. Spreader Bar</p>	<p style="text-align: center;">WP-SBAR-PL400</p>	<p style="text-align: center;">PL400E</p>		<p>It can be used on all PL lifts if don't mind the color difference.</p> <p>And will have different lifting range on all ProCare and HomeCare PL lifts</p>
 <p style="text-align: center;">Non-rotating 6 pt. Spreader Bar</p>	<p style="text-align: center;">WP-SBAR-PL600</p>	<p style="text-align: center;">PL500E PL600E</p>	<p style="text-align: center;">PL400E</p>	<p>It can be used on all PL lifts if don't mind the color difference.</p> <p>And will have different lifting range on all ProCare and HomeCare PL lifts</p>
 <p style="text-align: center;">2 Pt. Spreader Bar</p>	<p style="text-align: center;">DISCONTINUED</p>	<p style="text-align: center;">PL400E PL500E PL600E</p>		<p>Replace with WP-SBAR-PL400 or WP-SBAR-PL600</p>
 <p style="text-align: center;">Rotating 6 pt. Spreader Bar</p>	<p style="text-align: center;">DISCONTINUED</p>	<p style="text-align: center;">PL400E PL500E PL600E</p>		<p>Replace with WP-SBAR-PL400 or WP-SBAR-PL600</p>

Spreader Bar Compatibility Chart



Old Rotating 6 pt. Spreader Bar

DISCONTINUED

PL400E
PL500E
PL600E

Replace with WP-SBAR-PL400
or WP-SBAR-PL600